



# MRI probe materials



SWIFT, sw=64kHz,np/2=128, nv=4000, nt=16

**Not recommended as MRI probe materials:**

- 1(Acetal(POM)), 6(HIPS), 10(PVC/Acrylic), 11(Poliurethane, Thermost RIM),
- 13(PVC/Acrylic(white)), 14(PVC/Acrylic(blue)), 14pp(Polypropylene(white)),
- 15(RC 79D(yellow))

Acetal (polyoxymethylene copolymer)



PEEK (polyetheretherketone)

x16



Nylon 6/6 (polyimide)

x16




Ultem (polyetherimide)

x16




Polyphenylene Sulfide (polyphenylene sulfide)

x16



HIPS (high-impact polystyrene)



PTFE (polytetrafluoroethylene, mica filled)

x16


7



Ertalyte (polyethylene terephthalate)

x16

8



PVC/Acrylic (polyvinyl chloride/ polymethyl methacrylate)

10



Polyurethane (polyurethane, thermoset RIM)

11



Polyurethane (polyurethane, thermoset Casting)


x16

12




PVC/Acrylic (polyvinyl chloride/ polymethyl methacrylate)--white


13




PVC/Acrylic (polyvinyl chloride/ polymethyl methacrylate)--blue



14n--unkown




Polyurethane (RC-79D polyurethane)



Polyurethane (TD 277-02 polyurethane)



PVDF (Kunar)



## 1.0 BACKGROUND

High Field NMR requires that the background signal from the probe construction materials used in the device does not interfere with the signal from the sample. Subsequently it is necessary to perform an NMR analysis on all construction materials to be used within the probe.

## 2.0 PURPOSE AND SCOPE

In this experiment we will perform an NMR analysis on a variety of construction materials. We will also identify the materials suitable for High Field NMR and MRI applications. With these results we will produce a materials matrix that will contain information necessary for selecting the construction materials for NMR probes and RF coils for High Field MRI. The focus of the experiment is to verify which of the sample materials are suitable for the construction of a 3T  $^1\text{H}$  TEM head coil for commercial application.

## 3.0 REFERENCES

- E.E. Babcock, J.T. Vaughan, B. Lesan, R.L. Nunnally. Multinuclear NMR Investigations of Probe Construction Materials at 4.7T. Magn Reson Med 1990; 13: 498-503.

## 4.0 TEST SETUP

We selected 14 sample materials to be considered for NMR probe construction; they were supplied by independent manufacturers and fabricated into 3-in. lengths of 1-in. diameter rod. The testing will take place on the 9.4T/31cm Varian System at the University of Minnesota's Center for Magnetic Resonance Research (CMRR). The samples will be analyzed using five 1-in. diameter solenoid coils each tuned and matched to a frequency of interest:  $^1\text{H}$  (400.2 MHz),  $^{13}\text{C}$  (100.6 MHz),  $^{19}\text{F}$  (376.5 MHz),  $^{23}\text{Na}$  (105.9 MHz), and  $^{31}\text{P}$  (162.0 MHz).

## 5.0 TEST CASES

Each sample will be subjected to Free Induction Decay (FID) analysis covering the five established frequencies using the RF coils described above. The FID sequence and parameters will be determined by a qualified technician/scientist at the CMRR. Those results will be processed and recorded as NMR Spectra results for each sample material in the Test Results section (6.0).

After reviewing the data obtained from each sample we will determine which materials are suitable for probe construction, and the resulting definitions will be placed into the material matrix described above.

## 6.0 TEST RESULTS

### 6.1 Method / Test Parameters

The following sections to be completed during materials testing:

Sequence Description:

Sequence Parameters:

Signal Processing Description:

Misc. Description:

## 6.2 Sample Results

Table 1: Sum of FID

Materials	1H	13C	19F	23Na	31P
1: Acetal					
2: PEEK					
3: Nylon 6/6					
4: Ultem					
5: Polyphenylene Sulfide					
6: HIPS					
7: PTFE					
8: Ertalyte					
<del>10</del> 9: PVC/Acrylic: Royalite					
<del>11</del> 10: Polyurethane, RIM					
<del>12</del> 11: Polyurethane, Casting					
<del>13</del> 12: PVC/Acrylic, Kydex-100: White					
<del>14</del> 13: PVC/Acrylic, Kydex-100: Blue					
<del>15</del> 14: Polypropylene (white)					
15: RC 79D					
16: TD 277-02					
17:					

too large  
 ↓  
 15  
 16

too large  
 too large

Table 2: Presence (+) or Absence (-) of NMR Signal

Materials	1H	13C	19F	23Na	31P
1: Acetal					
2: PEEK					
3: Nylon 6/6					
4: Ultem					
5: Polyphenylene Sulfide					
6: HIPS					
7: PTFE					
8: Ertalyte					
<del>10</del> 9: PVC/Acrylic: Royalite					
<del>11</del> 10: Polyurethane, RIM					
<del>12</del> 11: Polyurethane, Casting					
<del>13</del> 12: PVC/Acrylic, Kydex-100: White					
<del>14</del> 13: PVC/Acrylic, Kydex-100: Blue					
<del>15</del> 14: Polypropylene					
15: RC 79D					
16: TD 277-02					
17:					

too large  
 ↓  
 15  
 16

Too large  
 Too large

Listed by id number, composition, and manufacturer are the NMR probe construction material candidates. Comment on the test results for each material's level of background signal found during the 9.4T experiment. Also recommend which sample materials could be used for the construction of a 3T <sup>1</sup>H single tuned TEM head coil for commercialization.

**Table 1: NMR probe constructions material candidates**

ID	Material	Composition	Manufacturer – Product Name
1	Acetal <i>1 sample</i>	Polyoxymethylene Copolmer (POM)	Quadrant EPP – Acetron GP, unfilled extruded
2	PEEK <i>1 sample</i>	Polyetheretherketone (PEEK)	Quadrant EPP – Ketron, PEEK 1000 unfilled extruded
3	Nylon 6/6 <i>1 sample</i>	Polyimide	Quadrant EPP – Nylon 101, type 66 nylon unfilled extruded
4	Ultem <i>1 sample</i>	Polyetherimide (PEI)	Quadrant EPP – Ultem 1000, unfilled extruded
5	Polyphenylene Sulfide <i>1 sample</i>	Polyphenylene Sulfide (PPS)	Quadrant EPP – Techron, unfilled extruded
6	HIPS <i>1 sample</i>	High-Impact PolyStyrene (HIPS)	TMX Ain Plastics
7	PTFE <i>1 sample</i>	Polytetrafluoroethylene (PTFE), Mica Filled	Quadrant EPP – Fluorosint 500, synthetic mica-filled, compression molded



8	Ertalyte 1 sample	Polyethylene Terephthalate (PET-P)	Quadrant EPP – Ertalyte, unreinforced, semi-crystalline thermoplastic polyester
9 10 on Bag	PVC/Acrylic 2 samples	Polyvinyl Chloride/Polymethyl Methacrylate	Uniroyal Technology Corp – Royalite R52, color: Ultra White
10 11 on Bag	Polyurethane 2 samples	Polyurethane, Thermoset RIM	Bayer AG Polyurethanes Div. - BAYDUR® 110 FR RIM Polyurethane, color: yellow
11 12 on Bag	Polyurethane 2 samples	Polyurethane, Thermoset Casting	Innovative Polymer – TD 277-02, color: white
12 13 on Bag	PVC/Acrylic 2 samples	Polyvinyl Chloride/Polymethyl Methacrylate	Kleerdex – KYDEX 100, color: White
13 14 on Bag	PVC/Acrylic 2 samples	Polyvinyl Chloride/Polymethyl Methacrylate	Kleerdex – KYDEX 100, color: Blue
14 15 No Bag	Polypropylene (white) 1 sample Labeled w/ Red Pen		Seelye – Eiler Plastics opaque white
15 16 Too large	RC 79D (yellow) 2 (Too large) samples		PDS Plastics

No sample #g →

Too large

Too large

16 16 Too large	TD 277-02 (white) 2 samples Too large		PDS Plastics
#55	PVDF 1 sample	KYNAR	white

<p>Experimenters:</p> <p>Test Date:</p>
---